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Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it

# MANUFACTURER INFORMATION ON REPROCESSING OF MEDICAL DEVICES ACCORDING TO EN ISO 17664:2017

**MANUFACTURER:** Larident Srl

**REF:** A11/A12/A13/A14

---

## WARNING NOTICE: NO SPECIAL NOTICE

---

**How to use it:** the medical device is used for pushing amalgama into the tooth. The medical device, amalgama syringe, is supplied in 4 different forms: Syringe with 60° curve (item code A11), Syringe with 10° curve (item code A12), Syringe with 90° curve (item code A13) and Syringe with 110° curve (item code A14).

**Storage and transport:** the product does not have special needs

**Cleaning:** take the device apart and clean the nozzle with a small brush

**Disinfection:** take the device apart and use suitable disinfectant solutions for disinfecting it

**Maintenance:** no special requirement

**Sterilization:** sterilize inside a sterilizing bag at 134° for 15 minutes, taking off the nozzle first. Autoclave cycles maximum: 50x

**Function Check and Test:** visual check for any damage or deformation of the product

**Storage:** the product does not have specific needs

The instructions listed above have been validated by the manufacturer of the medical devices as suitable for preparing the device for use.

It is the user's responsibility to ensure that the preparation process, the equipment the materials used and the staff obtain the desired results.

Any deviation from the recommended instructions must be carefully checked with regards to its efficiency and possible adverse consequences.

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**LARIDENT S.r.l.**
**PRODUCTEUR:** Larident Srl