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Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it

MANUFACTURER INFORMATION ON REPROCESSING OF MEDICAL DEVICES ACCORDING TO EN ISO 17664:2017

MANUFACTURER: Larident Srl

REF: A2/B A3/B

WARNING NOTICE: THE PRODUCT MUST BE COLD STERILIZED AT 121°

Limitation of reprocessing: 50x

How to use: medical device used to widen the cheek so as to have a better view of the operating field for intra-oral photography, for prophylaxis for orthodontic treatments

Storage and transport: no special requirement

Cleaning: use cleaning agent that takes into account the instructions for cleaning the product. After cleaning, remove cleanser thoroughly under running water. When disinfecting in thermodisinfectant do not exceed 93 °C

Disinfection: use disinfecting solutions that take into account the instructions for disinfecting the product, remove cleanser thoroughly under running water. When disinfecting in thermodisinfectant do not exceed 93 °C

Maintenance: no special requirement

Packaging: standardized packing material for sterilisation (recommendation)

Sterilisation: autoclavable vapor sterilisation at 121° for 15 minutes

Control / functional check: check the possible damages and possible distortion

Storage: no special requirement

The above-mentioned instructions have been validated as SUITABLE for preparation of medical device and its re-use by the medical device manufacturer.

It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing.

Normally, validation and routine controls are necessary.

Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

LARIDENT S.r.l.

Via Lamaineigra 12R (Ponte Rimassa)

16030 Tribogna (Genova) - ITALY

Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it