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Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it

MANUFACTURER INFORMATION ON REPROCESSING OF MEDICAL DEVICES ACCORDING TO EN ISO 17664:2017

MANUFACTURER: Larident Srl

REF: MT ASS.

WARNING NOTICE: THE PRODUCT MUST BE COLD STERILIZED

Limitation of reprocessing: single use

Storage and transport: no special requirement

Cleaning: cold sterilisation use cleaning agent that takes into account the instructions for cleaning the product. After cleaning, remove cleanser thoroughly under running water

Disinfection: use disinfecting solutions that take into account the instructions for disinfecting the product. remove cleanser thoroughly under running water

Maintenance: no special requirement

Packaging: the product must be packing in the package mini-grip

Sterilisation: cold sterilisation with disinfection liquids

Control / functional check: check the possible damages and possible distortion

Storage: no special requirement

The instructions listed above have been validated by the manufacturer of the medical devices as suitable for preparing the device for use.

It is the user's responsibility to ensure that the preparation process, the equipment the materials used and the staff obtain the desired results.

Any deviation from the recommended instructions must be carefully checked with regards to its efficiency and possible adverse consequences.

WATCH VIDEO

Scan or click

LARIDENT S.r.l.

Via Lamaineigra 12R (Ponte Rimassa)

16030 Tribogna (Genova) - ITALY

Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it