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Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it

# MANUFACTURER INFORMATION ON REUSING OF DM ACCORDING TO DIRECTIVE 93/42 CE

**MANUFACTURER:** Larident Srl

**REF:** R52

---

**WARNING NOTICE:** PROTECTIVE VISORS MADE OF POLYCARBONATE WITH AS INSERTED FRAME R52, ONLY AND EXCLUSIVELY PROTECT US FROM DROPS/SPRAYS OF WATER AND SALIVA

---

**Limitation of reusing:** the visor is reusable until its transparent effect is sufficient

**Assembling Instructions:**

- Remove all the protective films from the visor
- The visor must be positioned on the inside of the frame, with the nose opening towards the bottom
- Insert the visor into the metal tweezers of the frame
- Remember to adjust the visor according to the shape of your face so that they fit snugly on your cheekbones
- The visor is ready

**Storage and transport:** no special requirement

**Cleaning:** for the glasses shield use no alcohol cleaning, for the frame use cleaning agent that takes into account the instructions for cleaning the product. After cleaning remove cleanser thoroughly under running water

**Disinfection:** for the glasses shield use no alcohol disinfecting solutions, and the same for the frame. Remove cleanser thoroughly under running water

**Maintenance:** no special requirement

**Packaging:** the product must be packed in standard packaging

**Sterilisation:** frame not sterilizable

**Control / functional check:** check the possible damages and possible distortion

**Storage:** no special requirement

The above-mentioned instructions have been validated as SUITABLE for preparation of medical devices (DM) and its re-use by the medical device manufacturer.

It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing.

Normally, validation and routine controls are necessary.

Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

**LARIDENT S.r.l.**

Via Lamaneigra 12R (Ponte Rimassa)

16030 Tribogna (Genova) - ITALY

Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it