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
Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it

MANUFACTURER INFORMATION ON REUSING OF PPE ACCORDING TO DIN EN 166 2016/425

MANUFACTURER: Larident Srl

REF: R58

WARNING NOTICE: PROTECTIVE GLASSES MADE OF POLYCARBONATE WITH AS INSERTED FRAME R58, ONLY AND EXCLUSIVELY PROTECT AS FROM DROPS/SPRAYS OF WATER AND SALIVA

Limitation of reprocessing: the frame can be used until its transparency is adequate

Assembling Instructions:

- Remove all the protective films from the glasses
- The visor must be positioned on the inside of the frame, with the nose opening towards the bottom
- Widen the frame to insert the glasses' arms into two most external holes on each side
- Remember to adjust the glasses according to the shape of your face so that they fit snugly on your cheekbones
- The visor is ready. (see image 1)

Storage and transport: no special requirement

Cleaning: for the glasses shield use no alcohol cleaning, for the frame use cleaning agent that takes into account the instructions for cleaning the product. After cleaning, remove cleanser thoroughly under running water. When disinfecting in thermodisinfectant do not exceed 93 °C

Disinfection: for the glasses shield use no alcohol cleaning, and the same for the frame. Remove cleanser thoroughly under running water. When disinfecting in thermodisinfectant do not exceed 93 °C

Maintenance: no special requirement

Packaging: the product must be packed in the mini-grip package

Sterilisation: frame autoclavable vapor sterilisation at 121 °C

Control / functional check: check the possible damages and possible distortion

Storage: no special requirement


The above-mentioned instructions have been validated as SUITABLE for preparation of personal protection equipment (PPE) and its re-use by the medical device manufacturer.

It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing.

Normally, validation and routine check are necessary.

Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

image 1

LARIDENT S.r.l.

Via Lamaneigra 12R (Ponte Rimassa)

16030 Tribogna (Genova) - ITALY

Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it