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MANUFACTURER INFORMATION ON REPROCESSING OF (PPE) ACCORDING TO DIN EN 166 2016/425

MANUFACTURER: Larident Srl

REF: R62

WARNING NOTICE: PROTECTIVE VISORS MADE OF POLYCARBONATE WITH AS INSERTED FRAME R62, ONLY AND EXCLUSIVELY PROTECT FROM DROPS/SPRAYS OF WATER AND SALIVA

Limitation of reprocessing: the frame can be used until its transparency is adequate

Assembling Instructions:

- Remove all the protective films from the visor
- The visor must be positioned on the inside of the frame, with the nose opening towards the bottom
- Insert the visor into the metal tweezers of the frame
- Remember to adjust the visor according to the shape of your face so that they fit snugly on your cheekbones
- The visor is ready

Storage and transport: no special requirement

Cleaning: for the visor shield use no alcohol cleaning, for the frame use cleaning agent that takes into account the instructions for cleaning the product. After cleaning, remove cleanser thoroughly under running water. When disinfecting in thermodisinfector do not exceed 93 °C

Disinfection: for the visor shield use no alcohol disinfecting solutions, and the same for the frame. Remove cleanser thoroughly under running water. When disinfecting in thermodisinfector do not exceed 93 °C

Maintenance: no special requirement

Packaging: the product must be packed in the mini-grip package

Control / functional check: check the possible damages and possible distortion

Storage: no special requirement

The above-mentioned instructions have been validated as suitable for preparation of personal protection equipment (PPE) and its re-use by the medical device manufacturer.

It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing.

Normally, validation and routine checks are necessary.

Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

LARIDENT S.r.l.

Via Lamaineigra 12R (Ponte Rimassa)

16030 Tribogna (Genova) - ITALY

Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it