

## Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

**LARIDENT**


MADE IN ITALY  
SINCE 1985



### AVAILABLE COLORS



# TECHNICAL DATA SHEET

PRODUCT NAME	Bib holders
ITEM CODE	A33
MATERIAL/COMPONENTS	Santoprene (TPV)
THERMAL PROPERTIES	Autoclavable up to 121°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none"><li>• Package of 4 pcs. In carton box</li><li>• Available in 4 colors: white, fuchsia, yellow and green</li><li>• Plastic clip with spring</li></ul>
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines 
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

## Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

**LARIDENT**


MADE IN ITALY  
SINCE 1985



### GAMMA COLORI



# SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Catenelle in gomma termoplastica
NS CODICE	A33
MATERIALE	Santoprene (TPV)
PROPRIETÀ TERMICHE	Autoclavabile fino a 121°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none"><li>• Confezione da 4 pezzi in scatola di cartone</li><li>• Disponibile in 4 colori: bianco, fuxia, giallo e verde</li><li>• Clip di plastica con molla</li></ul>
PULIZIA	Utilizzare solamente detergenti adatti secondo le istruzioni del produttore
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente ad uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali 
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA