


Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

TECHNICAL DATA SHEET

| | |
|---|---|
| PRODUCT NAME | Cotton roll dispenser |
| ITEM CODE | A4 |
| MEASUREMENTS | cm 9,5 h10 |
| COLOUR | White, blue, yellow, grey |
| PACKAGING | 1 pc per pack. Packed individually in a cardboard box |
| MATERIAL | Coloured abs, nickel-plated iron weight, sealings in synthetic rubber |
| THERMAL PROPERTIES | Not sterilizable in an autoclave machine |
| MECHANICAL AND PHYSICAL CHARACTERISTICS | Contains 50 cotton rolls n.2. Can also be used with cotton rolls n.1 |
| CHEMICAL PROPERTIES | Resists all chemical disinfectants |
| SAFETY DATA SHEET | N/A |
| INSTRUCTIONS FOR USE | --- |
| PRECAUTIONS | Supplied in non-sterile conditions, before use clean the product |
| RISKS ANALYSIS IN REF. D.L. 46 24/2/97 | N/A |
| CLASSIFICATION REGULATION UE 2017/745 | It is not a medical device |

N/A = NOT APPLICABLE


N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

| | |
|---|---|
| PRODOTTO | Distributore di rulli a cassetto |
| NS CODICE | A4 |
| MISURE | cm 9,5 h10 |
| COLORE | Bianco, blu, giallo, grigio |
| CONFEZIONAMENTO | 1 pz per confezione. Imballo singolo in scatola di cartone |
| MATERIALE | Abs colorato, peso in ferro nichelato, guarnizioni in gomma sintetica |
| PROPRIETÀ TERMICHE | Non sterilizzabile in autoclave |
| PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE | Può contenere 50 rulli cotone del n.2. Si può usare anche con rulli del n.1 |
| PROPRIETÀ CHIMICHE | Resiste a tutti gli agenti chimici disinfettanti |
| SCHEDA DI SICUREZZA | N/A |
| ISTRUZIONI PER L'UTILIZZO | --- |
| PRECAUZIONI | Viene fornito in condizioni non sterili, prima dell'uso pulire l'articolo |
| ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97 | N/A |
| CLASSIFICAZIONE REGOLAMENTO UE 2017/745 | Non è un dispositivo medico |

N/A = NON APPLICABILE

N/N = NON NECESSARIA