

## Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it



### AVAILABLE COLORS



# TECHNICAL DATA SHEET

PRODUCT NAME	Cotton roll dispenser
ITEM CODE	A42
MATERIAL/COMPONENTS	Polyamide and polycarbonate
THERMAL PROPERTIES	Autoclavable up to 121°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none"><li>• Available in white or blue color</li><li>• The dispenser does not include cotton rolls</li><li>• Dimensions: 10,5x5 h 10 cm</li></ul>
CLEANING	Do not use alcohol for cleaning. Use specific cleaning products
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines 
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

## Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it



### GAMMA COLORI



# SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Portarulli a caduta individuale
NS CODICE	A42
MATERIALE	Poliammide e policarbonato
PROPRIETÀ TERMICHE	Autoclavabile a 121°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none"><li>• Disponibile nei colori bianco e blu</li><li>• La confezione non comprende rulli di cotone</li><li>• Dimensioni: 10,5x5 h 10 cm</li></ul>
PULIZIA	Non resiste all'alcol. Utilizzare agenti specifici
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente ad uso professionale per l'ergonomia dello studio
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali 
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA