

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

TECHNICAL DATA SHEET

PRODUCT NAME	Turboclean
ITEM CODE	A5
MATERIAL	Body: Polyamide Internal part: Acetal resin Discs: Foam-lined cardboard
THERMAL PROPERTIES	Autoclavable up to 134°C, absorbing discs not autoclavable
MECHANICAL AND PHYSICAL CHARACTERISTICS	• Measures: diameter: 12 cm, height: 12 cm. • Color: white • Container for tools dripping with 8 holes • Supplied with 5 absorbing discs
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE
N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEMA TECNICA DEL PRODOTTO

PRODOTTO	Turboclean
NS CODICE	A5
MATERIALE	Corpo: Poliammide Interno: Resina acetica Dischi: Cartone spugnato
PROPRIETÀ TERMICHE	Autoclavabile a 134°C, dischi assorbenti non autoclavabili
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	• Misure: diametro 12 cm, altezza 12 cm • Colore: bianco • Contenitore per sgocciolatura turbine, contrangoli, micromotori • 8 fori, fornito con 5 dischi assorbenti
PULIZIA	Utilizzare solamente detergenti adatti secondo le istruzioni del produttore
SCHEMA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative e alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A



N/A = NON APPLICABILE

N/N = NON NECESSARIA