

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



TECHNICAL DATA SHEET

PRODUCT NAME	Turbonet
ITEM CODE	A57
MATERIAL	Outer part: polypropylene loaded with calcium carbonate (marble powder) Brushes: nylon
THERMAL PROPERTIES	Material autoclavable at low speed (over 60° softens) it is necessary to remove the brush and suction pads (are not autoclavable)
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Measurements: cm 6,2x6,2 h 8• It is equipped with 8 brushes to clean the handles• It is equipped with 4 suction pads to keep it on the base stable• Equipped with a lid• Extra brushes are available cod. A57.1• For the dripping of the turbines the right container can be used cod. A5
CHEMICAL PROPERTIES	You can clean it with all the disinfectant products on the market
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	Autoclavable with attention to low speed after removing the brushes and suction pads
PRECAUTIONS	---
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A
CLASSIFICATION REGULATION UE 2017/745	Not a medical device

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Turbonet
NS CODICE	A57
MATERIALE	Corpo: polipropilene caricato con carbonato di calcio (polvere di marmo) Spazzole: nylon
PROPRIETÀ TERMICHE	Materiale autoclavabile a bassi regimi (oltre 60° si rammollisce) necessario smontare le spazzole e le ventose (non autoclavabili)
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Dimensioni: cm 6,2x6,2 h 8• Munito di 8 spazzole per la pulizia dei manipoli• Munito di 4 ventose sulla base per tenerlo stabile• Munito di coperchio• Sono disponibili le spazzole di ricambio cod. A57.1• Per la sgocciolatura delle turbine si può utilizzare l'apposito contenitore cod. A5
PROPRIETÀ CHIMICHE	Si può pulire con tutti i prodotti disinfettanti presenti sul mercato
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Autoclavare con attenzione ai bassi regimi dopo avere tolto le spazzole e le ventose
PRECAUZIONI	---
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico

N/A = NON APPLICABILE

N/N = NON NECESSARIA