


Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

TECHNICAL DATA SHEET

PRODUCT NAME	Hangable bib rolls dispenser
ITEM CODE	A64
MATERIAL/COMPONENTS	Polypropylene
THERMAL PROPERTIES	Not autoclavable
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Hangable bib rolls dispenser• Suitable for 11 cm diameter rolls max• Size: 28x13,5x12,5 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A




Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Contenitore per mantelline appendibile
NS CODICE	A64
MATERIALE	Polipropilene
PROPRIETÀ TERMICHE	Non autoclavabile
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Contenitore per mantelline appendibile• Adatto per un rotolo diametro max: 11 cm• Dimensioni: 28x13,5x12,5 cm
PULIZIA	Utilizzare solamente detergenti adatti secondo le istruzioni del produttore
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente ad uso professionale per l'ergonomia dello studio
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali 
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA