

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

TECHNICAL DATA SHEET

PRODUCT	Kidney shaped disposable cardboard dish
ITEM CODE	BIO-48
MATERIAL	100% recycled paper
THERMAL CHARACTERISTICS	N/A
MECHANICAL AND PHYSICAL CHARACTERISTICS	The device and packaging material are produced with non-toxic materials, without phthalates and latex. Light grey colour Internal measurements 250 x 140 x 45 mm, weight 25.6 g Approximately 700 ml capacity PAS29:1999 compliance - water test, held water for 24 hours (internal test)
INSTRUCTIONS FOR USE	Suitable for containing products to treat and/or clean the patient; suitable for any task that requires a disposable container. Product is not intended to be used as a container for hospital infected materials
WARNINGS	Keep out of direct sunlight and heat of any kind. Store at room temperature, protected from dust
DISPOSAL	---

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Bacinella reniforme monouso in cartone
NS CODICE	BIO-48
MATERIALE	Carta riciclata 100%
PROPRIETÀ TERMICHE	N/A
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	Il dispositivo e il materiale di confezionamento sono prodotti con materiali atossici esenti da ftalati e lattice. Colore grigio chiaro Misure interne 250 x 140 x 45 mm, peso 25,6 g Capacità 700 ml circa Conformità PAS29:1999 - prova con acqua tenuta 24 ore (test interno)
ISTRUZIONI PER L' UTILIZZO	Adatta a raccogliere materiale allo scopo di consentire la medicazione e/o pulizia del paziente; idonea per ogni operazione che richieda un contenitore monouso. Prodotto non destinato a essere utilizzato come contenitore di materiale ospedaliero infetto
PRECAUZIONI	Conservare al riparo da luce, fonti di calore, a temperatura ambiente e protetto da polvere
SMALTIMENTO	---

N/A = NON APPLICABILE

N/N = NON NECESSARIA