

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL DATA SHEET

PRODUCT	Plastic mixing pad made of PVC
OUR ITEM CODES	C10 - C11 - C12 - C13
MATERIAL	Non-toxic transparent PVC Non-slip base with adhesive paint Glue Latex free
THERMAL CHARACTERISTICS	Resistant up to 60° C
PHYSICAL AND MECHANICAL CHARACTERISTICS	The mixing pad is made up of 100 transparent PVC sheets each single sheet is 0,013 mm thick. The bottom of the pad is non-slip. C10 cm 5x5 pack of 6 pieces C11 cm 8x8 pack of 6 pieces C12 cm 14x8 pack of 3 pieces C14 cm 7x8 individually packaged Mixing pad with 50 sheets: C13 cm 14x20 pack of 3 pieces
CHEMICAL CHARACTERISTICS	We are not aware of intolerance of materials used in dental practices It cannot resist hydrocarbons Resistant to alcohol and oil
SAFETY DATA SHEET	N/A
WARNINGS	N/A
RISK ANALYSIS REF. D.L. 46 24/2/97	N/N

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Blocchetti per impasto in PVC
NS CODICE	C10 - C11 - C12 - C13 - C14
MATERIALE	PVC atossico trasparente Base con vernice antiscivolo Collante Latex free
PROPRIETÀ TERMICHE	Resistono fino a 60°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	Il blocchetto è composto da 100 fogli, il singolo foglio ha uno spessore di 0,013 mm. Il blocchetto ha fondo anti sdrucchiolo. C10 cm 5x5 conf. 6 pezzi C11 cm 8x8 conf. 6 pezzi C12 cm 14x8 conf. 3 pezzi C14 cm 7x8 impacchettato singolarmente Blocchetto composto da 50 fogli: C13 cm 14x20 conf. 3 pezzi
PROPRIETÀ CHIMICHE	Non siamo a conoscenza di intolleranza a materiali in uso nello studio dentistico. Non resiste a idrocarburi Resiste ad alcool ed olio
SCHEDA DI SICUREZZA	N/A
PRECAUZIONI	N/A
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/N

N/A = NON APPLICABILE

N/N = NON NECESSARIA