


Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



TECHNICAL SHEET

PRODUCT NAME	Disposable PET tray with raised edges and knurled base, matt white
ITEM CODE	C29
MATERIAL/COMPONENTS	PET
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Size 28x18 cm• Pack of 400 items divided into piles of 200 trays• With compartments for mixing liquids and divided inserts in the plastic for dental instruments
PRECAUTIONS/WARNINGS	Keep away from heat, smells and humidity. Do not store in the open air and in direct sunlight. Avoid heavy weights on the trays
CLASSIFICATION CEE 93/42, 47/2007	Not a medical device
Regulations/laws	Conforms to reach 1907/06 and additions Law 94/62/ce and additions

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Vassoio in PET monouso con bordino e fondo zigrinato, bianco opaco
NS CODICE	C29
MATERIALE	PET
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Dimensioni 28x18 cm• Confezione da 400 pezzi suddivisa in blocchi da 200• Con scomparti per pozzetti e porta strumenti
PRECAUZIONI	Tenere lontano da fonti di calore, odore e umidità non stoccare all'aperto e sotto il sole. Evitare la sovrapposizione di pesi eccessivi
CLASSIFICAZIONE CEE 93/42, 47/2007	Non è un dispositivo medico
NORME	Conforme a reach 1907/06 e succ.agg. Direttiva 94/62/ce e succ.agg.