

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



AVAILABLE COLORS



TECHNICAL DATA SHEET

PRODUCT NAME	Aligner box
ITEM CODE	G340 - G342 - G345
MATERIAL	Box: ABS Inner part: Polypropylene
THERMAL PROPERTIES	Cold disinfection
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• G340: without chewie clip• G342: with chewie clip• G345: with chewie clip and mirror• Available in 4 color combinations: blue/white, yellow/grey, black/grey and red/white• Dimensions: 10x7,6 h 2 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	N/N
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



GAMMA COLORI



Az/Bi G/Gr N/Gr Ro/Bi

G340

G342

G345

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Porta aligner
NS CODICE	G340 - G342 - G345
MATERIALE	Contenitore: ABS Interno: Polipropilene
PROPRIETÀ TERMICHE	Sterilizzare a freddo
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• G340: senza alloggiamento per chewie• G342: con alloggiamento per chewie• G345: con specchietto e alloggiamento per chewie• Disponibile in 4 combinazioni di colore: azzurro/bianco, giallo/grigio, nero/grigio e rosso/bianco• Misure: 10x7,6 h 2 cm
PULIZIA	Utilizzare solamente detergenti adatti secondo le istruzioni del produttore
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	N/N
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative e alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA