

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



G45 - G45 GLITTER




GLITTER



G46



G49



BITE BOX TECHNICAL PRODUCT SHEET

PRODUCT	Bite box with or without ventilation holes
ITEM CODE	G45 - G46 - G49
MATERIAL	Coloured non-toxic polypropylene
THERMAL PROPERTIES	Autoclavable up to 121°
MECHANICAL AND PHYSICAL PROPERTIES AND CHARACTERISTICS	The hinge is unbreakable Available in various colours Measurements: 7,8x7,3 h 2,5 cm (G45) - 7,8x7,7 h 3 cm (G46-G49) Sold in packs of ten items Packed in plastic bags
CHEMICAL PROPERTIES	Resists all chemical disinfectants
SAFETY SHEET	N/N
WARNINGS	Supplied non-sterile, clean before use
RISKS ANALYSIS IN REF. D. L.46 24/2/97, 37/2010	N/A
CLASSIFICATION Regulation (UE) 2017/745	It is not a medical device

N/A = NOT APPLICABLE
N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



G45 - G45 GLITTER




GLITTER



G46



G49



BITE BOX SCHEDA TECNICA

PRODOTTO	Bite box con o senza fori di ventilazione
CODICE ARTICOLO	G45 - G46 - G49
MATERIALE	Polipropilene atossico colorato
PROPRIETÀ TERMICHE	Autoclavabile fino a 121°
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	La cerniera è antirottura Disponibile in vari colori Misura: 7,8x7,3 h 2,5 cm (G45) - 7,8x7,7 h 3 cm (G46-G49) Venduti in confezione da 10 pezzi Confezionati in sacchetti di plastica
PROPRIETÀ CHIMICHE	Resiste bene ai disinfettanti chimici
SCHEDA DI SICUREZZA	N/N
PRECAUZIONI	Viene fornito in condizioni non sterili. Prima dell'uso pulire l'articolo
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97, 37/2010	N/A
CLASSIFICAZIONE Regolamento (UE) 2017/745	Non è un dispositivo medico

N/A = NON APPLICABILE

N/N = NON NECESSARIA