

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

TECHNICAL PRODUCT SHEET

PRODUCT	Box for cold disinfecting of dental instruments
ITEM CODE	G70
MATERIAL	polypropylene loaded with calcium carbonate (marble powder)
THERMAL CHARACTERISTICS	Autoclavable at low cycles (121°)
CHARACTERISTICS	Can contain about 3 litres of disinfectant. The plastic, holed basket that can be risen guarantees that the dental instruments dry off correctly. Measurements: cm 33x26 h14 Size of basket, internally: cm 27x17,50
CHEMICAL CHARACTERISTICS	Resistant to all chemical sterilizers
SAFETY DATA SHEET	N/A
PRECAUTIONS	None
ANALYSIS OF RISKS IN REF. D.L. 46 24/2/97	N/A
CLASSIFICATION CEE 93/42	It is not a medical device

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Box disinfezione a freddo
CODICE	G70
MATERIALE	Polipropilene caricato con carbonato di calcio (polvere di marmo)
PROPRIETÀ TERMICHE	Autoclavabile ai cicli più bassi (121°)
CARATTERISTICHE	Contiene circa 3 litri di disinfettante Il cestello sollevandosi garantisce l'asciugatura degli strumenti Misure: cm 33x26 h14 Misure utili del cestello interno: cm 27x17,50
PROPRIETÀ CHIMICHE	Resistente a tutti gli sterilizzanti chimici
SCHEDA DI SICUREZZA	N/A
PRECAUZIONI	Nessuna
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A
CLASSIFICAZIONE CEE 93/42	Non è un dispositivo medico

N/A = NON APPLICABILE

N/N = NON NECESSARIA