

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



AVAILABLE COLORS



TECHNICAL DATA SHEET

PRODUCT NAME	Round burs holder
ITEM CODE	L47
MATERIAL/COMPONENTS	Base: Polyamide Lid: Polycarbonate
THERMAL PROPERTIES	Autoclavable up to 121°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Round burs holder completely sterilizable at 121°C, with anti-shred base• Measures: ø 8 h 9,5 cm• Available in white or grey• Available in the variants: FG, RA and Mixed
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



GAMMA COLORI



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Portafrese rotondo
NS CODICE	L47
MATERIALE	Base: Poliammide Coperchio: Policarbonato
PROPRIETÀ TERMICHE	Autoclavabile a 121°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Portafrese rotondo completamente sterilizzabile a 121°C, con base antistrucchio• Misure: ø 8 h 9,5 cm• Disponibile bianco o grigio• Disponibile nelle varianti: FG, RA e Misto
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA