

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL DATA SHEET

PRODUCT NAME	Burs holder
ITEM CODE	LS12 - LS18 - LS24 - LS36 - LS48 - LS72
MATERIAL/COMPONENTS	Base: Polypropilene with marmor powder Lid: Polystirene
THERMAL PROPERTIES	Autoclavable base up to 134°C. Non autoclavable lid
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Bursholder available in the formats: FG, RA, Misto or HP• 12, 18, 24, 36, 48 e 72 holes• Measures: 12 to 36 holes: 6,5 x 5 h 3 cm with HP lid h 5,5 cm• Measures: 48 and 72 holes: 8,5 x 6,2 h 3 cm with HP lid h 5,5 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE


N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Portafrese
NS CODICE	LS12 - LS18 - LS24 - LS36 - LS48 - LS72
MATERIALE	Base: Polipropilene con polvere di marmo Coperchio: Polistirolo
PROPRIETÀ TERMICHE	Base autoclavabile a 134°C. Coperchio non autoclavabile
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none"> • Portafrese disponibile nei formati: FG, RA, Misto o HP • 12, 18, 24, 36, 48 e 72 fori • Misure: 12 a 36 fori: 6,5 x 5 h 3 cm con coperchio HP: h 5,5 cm • Misure: 48 e 72 fori: 8,5 x 6,2 h 3 cm con coperchio HP: h 5,5 cm
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA