

**Larident Srl**

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

**LARIDENT**MADE IN ITALY  
SINCE 1985

## TECHNICAL DATA SHEET

PRODUCT NAME	Endobox
ITEM CODE	Q15
MATERIAL/COMPONENTS	Aluminium
THERMAL PROPERTIES	Autoclavable up to 134°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none"><li>• Mini endobox mixed holes, red</li><li>• Measures: 10x4,2 h 5 cm</li></ul>
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

**Larident Srl**

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

**LARIDENT**MADE IN ITALY  
SINCE 1985

## SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Endobox
NS CODICE	Q15
MATERIALE	Alluminio
PROPRIETÀ TERMICHE	Autoclavabile a 134°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none"><li>• Mini endobox fori misti, rosso</li><li>• Misure: 10x4,2 h 5 cm</li></ul>
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA