

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL DATA SHEET

PRODUCT NAME	Tray rack
ITEM CODE	Q19 - Q20 - Q24
MATERIAL/COMPONENTS	Aluminium
THERMAL PROPERTIES	Autoclavable up to 134°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Tray rack available in various sizes and models:<ul style="list-style-type: none">- Q19: for 11 instruments (for standard tray 28x18 cm)- Q20: for 17 instruments (for standard tray 28x18 cm)- Q24: for 8 instruments (for small tray 14x18 cm)
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Rastrelliera per vassoi
NS CODICE	Q19 - Q20 - Q24
MATERIALE	Alluminio
PROPRIETÀ TERMICHE	Autoclavabile fino a 134°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Rastrelliera per vassoi disponibile in varie dimensioni e modelli:<ul style="list-style-type: none">- Q19: per 11 strumenti (per vassoio standard 28x18 cm)- Q20: per 17 strumenti (per vassoio standard 28x18 cm)- Q24: per 8 strumenti (per vassoio piccolo 14x18 cm)
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA