

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595


info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



AVAILABLE COLORS



TECHNICAL DATA SHEET

PRODUCT NAME	Perforated aluminium tray
ITEM CODE	Q34
MATERIAL/COMPONENTS	Anodized aluminium
THERMAL PROPERTIES	Autoclavable up to 134°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">Perforated anodized aluminium tray available in the colors: aluminium, blue, yellow, red and greenMeasures: 18x14 h 1,5 cm
CLEANING	Use a cleaning agent suitable for the material, not suitable to use with a thermodisinfectant as the color can fade
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



GAMMA COLORI



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Vassoio in alluminio forato
NS CODICE	Q34
MATERIALE	Alluminio anodizzato
PROPRIETÀ TERMICHE	Autoclavabile a 134°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Vassoio in alluminio anodizzato forato disponibile nei colori: alluminio, blu, giallo, rosso e verde• Misure: 18x14 h 1,5 cm
PULIZIA	Utilizzare solamente detergenti adatti al materiale, non adatto al termodisinfettore in quanto potrebbe scolorire
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA