

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL DATA SHEET

PRODUCT NAME	Container for endodontic instruments
ITEM CODE	Q62
MATERIAL/COMPONENTS	Polypropylene with marmor powder
THERMAL PROPERTIES	Autoclavable up to 121°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Endo Clean Test: round container for 16 endodontic instruments, with graduated length gauge and sponge for cleaning• Size: ø7,4 h 6,8 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Contenitore per strumenti canalari
NS CODICE	Q62
MATERIALE	Polipropilene caricato marmo
PROPRIETÀ TERMICHE	Autoclavabile fino a 121°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Endo Clean Test: contenitore rotondo per 16 strumenti canalari, con misuratore a scalare per lunghezza, spugnetta per pulizia• Misure: $\varnothing 7,4$ h 6,8 cm
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA