

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



TECHNICAL DATA SHEET

PRODUCT NAME	Dappen pots
ITEM CODE	R16
MATERIAL/COMPONENTS	Mixing pots: Santoprene Lid: Polystyrene
THERMAL PROPERTIES	Autoclavable up to 121°C non autoclavable lid
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Mixing pots with orange UV protective lid• 5 bases with 4 cavities + 1 lid• Measures: base 5x5 pots ø 1,9 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Pozzetti per miscelazione
NS CODICE	R16
MATERIALE	Pozzetti: Santoprene Coperchio: Polistirolo
PROPRIETÀ TERMICHE	Autoclavabile a 121°C coperchio non autoclavabile
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Pozzetti per miscelazione con coperchio arancione protettivo da raggi UV• 5 basi da 4 cavità + 1 coperchio• Misure: base 5x5 pozzetto ø 1,9 cm
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA