


Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



TECHNICAL DATA SHEET

PRODUCT NAME	Spatula for cement and acrylic
ITEM CODE	R23 - R24
MATERIAL/COMPONENTS	Polyamide
THERMAL PROPERTIES	Autoclavable up to 121°C
MECHANICAL AND PHYSICAL CHARACTERISTICS	<ul style="list-style-type: none">• Polyamide spatula for cement and acrylic- R23: package 10 pcs- R24: package 100 pcs• Measures: 18,5x1,2 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE


N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Spatola per cementi e resine
NS CODICE	R23 - R24
MATERIALE	Poliamide
PROPRIETÀ TERMICHE	Autoclavabile a 121°C
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	<ul style="list-style-type: none">• Spatola in poliamide per cementi resine acriliche e paste- R23: conf. 10 pz- R24: conf. 100 pz• Misure: 18,5x1,2 cm
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA