

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL SHEET

PRODUCT NAME	PROTECTIVE VISOR IN PET, WITH SPONGE AND ELASTIC
ITEM CODE	R4
COMPONENTS	PROTECTIVE VISOR: PET VISOR THICKNESS: 0.300 mm SPONGE: HCL30 IP/PAT
THERMAL CHARACTERISTICS	THE PRODUCT IS SINGLE USE AND NOT STERILISABLE
PHYSICAL CHARACTERISTICS, MECHANICAL CHARACTERISTICS AND PACKAGING	PACKAGING IN CARTON BOX: 12 PCS THE VISOR COVERS NOSE AND MOUTH THE PROTECTIVE VISOR HAS A DOUBLE PROTECTIVE COATING THE VISOR IS VERY TRANSPARENT
CHEMICAL CHARACTERISTIC	ONLY DISINFECT WITH NON-ALCOHOLIC SOLUTIONS.
SAFETY DATA SHEET	N/A
PRECAUTION	WARNING: ONLY USE AS FACE SHIELD AGAINST SPLASHING LIQUID OF WATER SALIVA.
INSTRUCTIONS FOR USE	ONLY TO BE USED BY MEDICAL SPECIALISTS OR HIGHLY SPECIALISED STAFF. PRODUCED NON-STERILE.
ANALYSIS RISKS REF.D.L.46 24/2/97	ONLY AVAILABLE IN ITALIAN
CLASS CEE EN 166	IS AN INDIVIDUAL PROTECTION DEVICE (DPI)

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	VISIERA PROTETTIVA IN PET, CON SPUGNA ED ELASTICO
NS CODICE	R4
MATERIALE	VISIERA: PET SPESSORE VISIERA : 0.300 mm SPUGNA: HCL30 IP/PAT
PROPRIETA' MECCANICHE, FISICHE, CARATTERISTICHE E CONFEZIONAMENTO	CONFEZIONE: 12 PEZZI. LA VISIERA CHE COPRE NASO E BOCCA, HA UN DOPPIO RIVESTIMENTO PROTETTIVO ED E' DOTATA DI UNA ECCELLENTE TRASPARENZA
PROPRIETA' CHIMICHE	DISINFETTARE SOLAMENTE LA VISIERA CON SOLUZIONI NON ALCOLICHE
SCHEDA DI SICUREZZA	N/A
ISTRUZIONI PER L' UTILIZZO	DESTINATI AD ESSERE UTILIZZATI ESCLUSIVAMENTE DAL MEDICO SPECIALISTA O DA PERSONALE ALTAMENTE SPECIALIZZATO. VIENE FORNITO IN CONDIZIONE "NON STERILE"
PRECAUZIONI	ATTENZIONE: LA VISIERA PROTETTIVA PROTEGGE SOLO ED ESCLUSIVAMENTE DA GOCCE E SPRUZZI D'ACQUA O SALIVA (DROPLET)
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97, 37/2010	SI
CLASSIFICAZIONE secondo la norma UNI EN 166:2004 e il regolamento 2016/425	E' UN DISPOSITIVO DI PROTEZIONE INDIVIDUALE (DPI) DI CAT II (DPI)

N/A = NON APPLICABILE

N/N = NON NECESSARIA