

## Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

**LARIDENT**

MADE IN ITALY  
SINCE 1985



# TECHNICAL SHEET

PRODUCT NAME	PROTECTIVE EYESHIELDS WITH AUTOCLAVABLE FRAME
ITEM CODE	R58
COMPONENTS	PROTECTIVE VISOR: PET VISOR THICKNESS: 0.300 mm FRAME: MAKROLON 2458 (POLYCARBONATE)
THERMAL CHARACTERISTICS	FRAME AUTOCLAVABLE UP TO 121°C.
PHYSICAL CHARACTERISTICS, MECHANICAL CHARACTERISTICS AND PACKAGING	THE FRAME IS AVAILABLE IN 5 COLOURS. PACKAGING: 1 FRAME + 30 DISPOSABLE LENS WITH MOVABLE LENS. THE VISOR HAS A DOUBLE PROTECTIVE FILM.
CHEMICAL CHARACTERISTICS	IT RESISTS TO THE CHEMICAL DISINFECTANTS. ONLY DISINFECT WITH NON-ALCOHOLIC SOLUTIONS.
SAFETY DATA SHEET	N/A
PRECAUTION	WARNING: ONLY USE AS FACE SHIELD AGAINST SPLASHING LIQUID OF WATER SALIVA.
INSTRUCTIONS FOR USE	ONLY TO BE USED BY MEDICAL SPECIALISTS OR HIGHLY SPECIALISED STAFF. SUPPLIED "NON-STERILE".
ANALYSIS RISKS REF.D.L.46 24/2/97	ONLY AVAILABLE IN ITALIAN.
CLAS. CEE EN 166	IS AN INDIVIDUAL PROTECTION DEVICE (DPI)

## Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

**LARIDENT**

MADE IN ITALY  
SINCE 1985



# SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	OCCHIALI PROTETTIVI IN POLICARBONATO CON MONTATURA AD INCASTRO AUTOCLAVABILE.
NS CODICE	R58
MATERIALE	VISIERA: PET SPESSORE VISIERA : 0.300 mm MONTATURA: MAKROLON 2458 (POLICARBONATO)
PROPRIETÀ TERMICHE	LA MONTATURA E' STERILIZZABILE A 121°C.
PROPRIETÀ MECCANICHE, FISICHE, CARATTERISTICHE E CONFEZIONAMENTO	LA MONTATURA E' DISPONIBILE IN 5 COLORI. CONFEZIONE: 1 MONTATURA + 30 OCCHIALI MONOUSO CON LENTI ORIENTABILI. LA VISIERA HA UN DOPPIO RIVESTIMENTO PROTETTIVO.
PROPRIETÀ CHIMICHE	DISINFETTARE SOLAMENTE LA VISIERA CON SOLUZIONI NON ALCOLICHE
SCHEDA DI SICUREZZA	N/A
ISTRUZIONI PER L'UTILIZZO	DESTINATO AD ESSERE UTILIZZATO ESCLUSIVAMENTE DAL MEDICO SPECIALISTA O DA PERSONALE ALTAMENTE SPECIALIZZATO. VIENE FORNITO IN CONDIZIONE "NON STERILE".
PRECAUZIONI	ATTENZIONE: GLI OCCHIALI PROTETTIVI PROTEGGONO SOLO ED ESCLUSIVAMENTE DA GOCCE E SPRUZZI D'ACQUA O SALIVA (DROPLET).
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97, 37/2010	SI
CLASSIFICAZIONE secondo la norma UNI EN 166:2004 e il regolamento 2016/425	E' UN DISPOSITIVO DI PROTEZIONE INDIVIDUALE (DPI)

N/A = NON APPLICABILE

N/N = NON NECESSARIA