

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL SHEET

PRODUCT NAME	ACRILIC PROTECTIVE VISOR WITH AUTOCLAVABLE 121° INTEGRATED CLIP FRAME
CODE	R60
COMPONENT	PROTECTIVE VISOR: LUMIRROR VISOR THICKNESS: 0.175 mm FRAME: MAKROLON 2458 (POLYCARBONATE)
THERMAL CHARACTERISTIC	THE VISOR IS REUSABLE AND NOT STERILISABLE
PHYSICAL CHARACTERISTICS, MECHANICAL CHARACTERISTICS AND PACKING	THE SHIELD R60 IS EXTRA-THIN 0.175 mm THE SHIELD IS NOT SCRATCHABLE AND ANTI-FOG THE R60 PACKED INCLUDING: 6 VISOR AND 1 FRAME SHIELDS
CHEMICAL CHARACTERISTIC	IT RESISTS TO THE CHEMICAL DISINFECTANTS
SAFETY DATA SHEET	N/A
PRECAUTION	WARNING: USE ONLY AS FACE SHIELD AGAINST SPLASHING LIQUID OF WATER SALIVA
INSTRUCTIONS FOR USE	SUPPLIED IN "NON-STERILE" CONDITION. AUTOCLAVABLE. THE PRODUCT IS SINGLE USE. AVAILABLE INSTRUCTION FOR FITTING INSIDE CARTON BOX.
ANALYSIS OF THE RISKS REF.D.L.46 24/2/97	AVAILABLE ONLY IN ITALIAN
CLAS. CEE EN 166	IS AN INDIVIDUAL PROTECTION DEVICE (DPI)

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	VISIERA PROTETTIVA IN MATERIALE ACRILICO CON MONTATURA AD INCASTRO AUTOCLAVABILE A 121°
NS CODICE	R60
MATERIALE	VISIERA: LUMIRROR SPESSORE VISIERA: 0.175 mm FRAME: MAKROLON 2458 (POLICARBONATO)
PROPRIETÀ TERMICHE	LO SCHERMO PROTETTIVO È RIUTILIZZABILE E NON STERILIZZABILE
PROPRIETÀ MECCANICHE, FISICHE, CARATTERISTICHE, CONFEZIONAMENTO	CONFEZIONE: 1 MONTATURA + 6 VISIERE LA VISIERA R60 EXTRA SOTTILE 0.175 mm, ANTISTATICA, ANTIGRAFFIO, ANTI APPANNAMENTO
PROPRIETÀ CHIMICHE	RESISTE A TUTTI GLI AGENTI CHIMICI DISINFETTANTI
SCHEDA DI SICUREZZA	N/A
PRECAUZIONI	ATTENZIONE: LA VISIERA PROTETTIVA PROTEGGE SOLO ED ESCLUSIVAMENTE DA GOCCE E SPRUZZI D'ACQUA O SALIVA
ISTRUZIONI PER L'UTILIZZO	DESTINATA AD ESSERE UTILIZZATA ESCLUSIVAMENTE DAL MEDICO SPECIALISTA O DA PERSONALE ALTAMENTE SPECIALIZZATO. VIENE FORNITA IN CONDIZIONE "NON STERILE" ED IL PRODOTTO È MONOUSO
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97, 37/2010	SI
CLASSIFICAZIONE CEE EN 166	E' UN DISPOSITIVO DI PROTEZIONE INDIVIDUALE (DPI)

N/A = NON APPLICABILE

N/N = NON NECESSARIA