

Larident Srl


Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



AVAILABLE COLORS



TECHNICAL DATA SHEET

PRODUCT NAME	Dappen
ITEM CODE	V2
MATERIAL/COMPONENTS	Glass
THERMAL PROPERTIES	Cold disinfection. If autoclaving is necessary, proceed with caution and allow to cool slowly
MECHANICAL AND PHYSICAL CHARACTERISTICS	• Glass dappen • Available in colours: amber, blue, green and transparent • Measures: ø 2 h 3 cm
CLEANING	Use a cleaning agent suitable for the material
SAFETY DATA SHEET	N/N
INSTRUCTIONS FOR USE	This product is for professional use only
PRECAUTIONS/WARNINGS	Non-sterile
DISPOSAL	Disposal of the product must comply with local regulations and environmental guidelines
CLASSIFICATION REGULATION UE 2017/745	Not a medical device
RISKS ANALYSIS IN REF. D.L. 46 24/2/97	N/A

N/A = NOT APPLICABLE


N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

**GAMMA COLORI**

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Dappen
NS CODICE	V2
MATERIALE	Vetro
PROPRIETÀ TERMICHE	Disinfezione a freddo. Se necessaria autoclavazione prestare attenzione e raffreddare lentamente
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	• Dappen in vetro • Disponibile nei colori: ambra, blu, verde e trasparente • Misure: ø 2 h 3 cm
PULIZIA	Utilizzare solamente detergenti adatti al materiale
SCHEDA DI SICUREZZA	N/N
ISTRUZIONI PER L'UTILIZZO	Destinato esclusivamente a uso professionale
PRECAUZIONI	Non sterile
SMALTIMENTO	Lo smaltimento del prodotto deve essere conforme alle normative ed alle linee guida ambientali locali
CLASSIFICAZIONE REGOLAMENTO UE 2017/745	Non è un dispositivo medico
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A

N/A = NON APPLICABILE

N/N = NON NECESSARIA